

COVID-19 Antigen Rapid Test Cassette (swab) Clinical Sensitivity and Specificity Study Report

1. Objective

The HEO[®] COVID-19 Antigen Rapid Test Cassette (swab) manufactured by Hangzhou HEO Technology Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the HEO COVID-19 Antigen Rapid Test Cassette (swab) and the comparator RT-PCR assay.

2. Method

A study of 628 direct nasopharyngeal swabs specimens was performed. The specimens were prospectively collected from patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care at 5 locations and tested at a single central laboratory.

Two nasopharyngeal swab specimen were collected from individual symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19. At all locations, one nasopharyngeal swab specimen was tested directly with the COVID-19 Antigen Rapid Test Cassette (swab) according to product instructions for use, and another nasopharyngeal swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. Nasopharyngeal swabs specimens were tested by operators who were blinded to the RT-PCR test result.

The positive percent agreement (PPA) was calculated as $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$. The negative percent agreement (NPA) was calculated as $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$. The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd., is a real-time reverse transcription polymerase chain reaction (RT-PCR) test. This product has got CE, NMPA certifications and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37 and the Ct value of human housekeeping gene β -Actin is not higher than 35.

4. Enrollment criteria (inclusion/exclusion criteria)

4.1 Inclusion criteria

- Patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care.
- Symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19.

4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

5. Result

The results are summarized in the following table.

The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The

sensitivity was calculated for the different Ct value range (Ct value \leq 37)

COVID-19 Antigen		RT-PCR (Ct value \leq 37)		Total
		Positive	Negative	
HEO[®]	Positive	172	0	172
	Negative	3	453	456
Total		175	453	628

PPA (Ct \leq 37):98.28% (172/175), (95%CI: 95.08%~99.41%)
 NPA: 100% (453/453), (95%CI: 99.15%~100%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

6. Conclusion

Taken together, the HEO[®] COVID-19 Antigen Rapid Test Cassette (Lab) had a positive percent agreement (sensitivity) of 98.28% (95% CI: 95.08%~99.41%) with specimens of a Ct count \leq 37, and negative percent agreement (specificity) of 100% (95% CI: 99.15%~100%).

